



Agilia® SP TIVA (TCI)

Syringe infusion devices for intravenous anaesthesia





Agilia SP TIVA

Where IV meets Anaesthesia

Agilia SP TIVA is the syringe device dedicated to Anaesthesia

Agilia SP TIVA with the following pharmacokinetic models, all in plasma as well as effect target:

- Marsh & Schnider for Propofol administration in adults
- **Kataria & Paedfusor** for administration in paediatrics
- Minto, Gepts & Scott for Remifentanil, Sufentanil and Alfentanil administration in adults

The Agilia SP TIVA provides pharmacokinetic TCI models covering the widest range of patient population, including children



Agilia Connect is as intuitive as ever and now comes with embedded WiFi

Latest standards enables cyber secure communication



Intuitive setting


- The infusion setting is guided and each step clearly identified. Modes mL/h, Dose and TCI at your disposal
- The displayed instructions are simple to follow. Intuitive settings increase user confidence





Clear monitoring

- Monitoring is simple thanks to easily accessible numerical and graphical information
- The syringe is always visible and the flow indicator lights enable quick and safe monitoring



Quick history access



- TCI Infusion history, both in dose and flow rate
- Key infusion parameters
- Target and concentration history
- · Wake up time
- In TIVA mode, dose rate or mL/h, a wide flexibility to cover your whole library from your TIVA device
- Assisted with Vigilant Master Med, Agilia SP TIVA brings you a higher level of safety in Total Intravenous Anaesthesia



